

Average dimensional change measured in mm at room temperature (19C) inside calibrated glass tubes.

| Nominal Outer Diameter x Length | Expanded Diameter | Expanded Length | Nominal Outer Diameter x Length | Expanded Diameter | Expanded Length |
|---------------------------------|-------------------|-----------------|---------------------------------|-------------------|-----------------|
| 8.5/10 x 60 | 5 | 81 | 10/12 x 90 | 7 | 123 |
| | 7 | 77 | | 8 | 118 |
| | 8 | 72 | | 9 | 113 |
| 8.5/10 x 70 | 5 | 95 | 12/14 x 60 | 7 | 86 |
| | 7 | 89 | | 9 | 81 |
| | 8 | 84 | | 11 | 75 |
| 8.5/10 x 80 | 5 | 108 | 12/14 x 70 | 7 | 100 |
| | 7 | 102 | | 9 | 94 |
| | 8 | 96 | | 11 | 87 |
| 10/12 x 60 | 7 | 83 | 12/14 x 80 | 7 | 114 |
| | 8 | 79 | | 9 | 108 |
| | 9 | 76 | | 11 | 100 |
| 10/12 x 70 | 7 | 96 | 12/14 x 90 | 7 | 128 |
| | 8 | 92 | | 9 | 121 |
| | 9 | 88 | | 11 | 112 |
| 10/12 x 80 | 7 | 110 | | | |
| | 8 | 105 | | | |
| | 9 | 101 | | | |