

Aesculap Bipolar Cables



Aesculap Neurosurgery

Bipolar Cables

Aesculap offers a wide range of different bipolar instruments like YASARGIL bipolar forceps, CASPAR or SINTRAM NON-STICK bipolar forceps which are tailored for neurosurgical treatments. Additionally provided is a bipolar unit GN060, an irrigation unit GN090 plus all accessories like bipolar cables for a complete set for the operating room.


To use our bipolar instruments in combination with other units, we can offer several bipolar cables which have either a 22 mm pin distance – suitable e.g. for Dolly or EMC units, or a 38 mm pin distance – suitable e. g. for Codman units. Furthermore, round pin connector or a special connector for Codman units complete the program.

All the bipolar cables fulfill the new standard of Edition 4 of IEC-60601-2-2.


Connector to the electrosurgical unit

Equipment side


PIN Distance 22 mm
EMC, Dolley



PIN Distance 28.6 mm
Aesculap, BOWA, EMC,
Lamidey, Söring, Valleylab,
ERBE-International



PIN Distance 38 mm
Codman CMC II, CMC III



5 mm round connector, 2 mm PIN
Aesculap, ACMI, Berchtold, Martin, Wolf



8 mm round connector, 4 mm PIN
ERBE, Storz



Connector to the instrument

Instrument side

 Aesculap tab connector

 Aesculap Euro tab connector

 Aesculap round pin connector



Equipment side	Instrument side	Article-No.	Pin Distance	Description
		•• GK331	28.6 mm	4 m length, suitable for units like: Aesculap, BOWA, EMC, Lamidey, Söring, Valleylab
		•• GK332	38 mm	4 m length, suitable for units like: Codman CMC II, CMC III
		GK342	ø 8 mm/ ø 4 mm	4 m length, suitable for units like: ERBE, Storz
		GN140	ø 5 mm/ ø 2 mm	Round pin, 4 m length, suitable for units like: Aesculap, ACMI, Berchtold, Martin, Wolf
		•• GN141	ø 5 mm/ ø 2 mm	Round pin, 4 m length, suitable for units like: Aesculap, ACMI, Berchtold, Martin, Wolf
		GK195	28.6 mm	4 m length, suitable for units like: Aesculap, BOWA, EMC, Lamidey, Söring, Valleylab
		GN145	22 mm	5 m length, suitable for units like: Dolly, EMC, ERBE-International
		GN146	38 mm	5 m length, suitable for units like: Codman CMC II, CMC III

Equipment side	Instrument side	Article-No.	Pin Distance	Description
		 GK194	28.6 mm	4 m length, suitable for units like: Aesculap, BOWA, EMC, Lamidey, Söring, Valleylab
		 GN076	ø 5 mm/ ø 2 mm	4 m length, suitable for units like: Aesculap, ACMI, Berchtold, Martin, Wolf
		 GK185	ø 8 mm/ ø 4 mm	4 m length, suitable for units like: ERBE, Storz
		 GK217		4 m length, suitable for units like: Codman

Equipment side	Instrument side	Article-No.	Pin Distance	Description
		 GK280	22 mm	4 m length, suitable for units like: Dolly, EMC, ERBE-International
		 GK281	28.6 mm	4 m length, suitable for units like: Aesculap, BOWA, EMC, Lamidey, Söring, Valleylab
		 GK282	38 mm	4 m length, suitable for units like: Codman CMC II, CMC III
		 GK290SU 	22 mm	3.7 m length, suitable for units like: Dolly, EMC, ERBE-International
		 GK291SU 	28.6 mm	3.7 m length, suitable for units like: Aesculap, BOWA, EMC, Lamidey, Söring, Valleylab
		 GK292SU 	38 mm	3.7 m length, suitable for units like: Codman CMC II, CMC III
		 GK328	UHF	4 m length, suitable for units like: Aesculap GK050, GK051
		 GK329	UHF	4 m length, suitable for units like: Aesculap GK050, GK051

