

PRODUCT DATA SHEET

EXADROP(4061209)

Nombre del Atributo	Valor de Atributo
Marca	Exadrop®
Product name (HC)	Exadrop®
Tipo de producto	Administration set
Application	Infusion
Type of administration set	Gravity
Container connector	IV standard spike
Longitud	150.00 cm
Inner-ø tube	3.00 mm
Outer-ø tube	4.10 mm
Type of filter	Fluid filter
Position of filter	Drip chamber
Pore size (Fluid filter)	15 µm
drops/ml	20
Air ventilation	vented
Flow control device	Regulator, Slide clamp
with AirStop	No
with PrimeStop	No

Nombre del Atributo	Valor de Atributo
Number of back check valves	-
Number of needle-free valves	-
Injection access	-
Connection	Luer-Lock/male
Light-protected	No
Color of light protection	-
Pump compatibility (Type)	-
Material	PVC
Tipo de packaging	Peel pack
Producto (sin empacar) contie	Si
producto contiene LATEX	No
Type of Sterilizat. (Current)	No
Producto (sin empacar) contie	No
producto contiene LATEX	No
Type of Sterilizat. (Current)	No
producto contiene PVC	Yes
producto contiene LATEX	No
producto contiene DEHP	No
producto contiene níquel	No
Type of Sterilizat. (Current)	Estéril
con rosca/suave	EO