

Aesculap[®] Trocar Solution

Reusable and Hybrid Systems



Aesculap Endoscopic Technology

Aesculap® Trocar Solution



— Hybrid Solution —

Combines the advantages of reusable trocars with the high functionality of single-use components.





- low instrument friction
- firm hold in the abdominal wall
- easy to clean, as the system can be fully disassembled

Simple

- excellent tactile feedback during insertion into the abdomen
- low space requirement through slender trocar design
- exceptional ergonomic properties and good weight distribution

Exceptional

Modular

- individual combination of the system based on a reusable trocar sleeve
- choice between a single-use universal seal and a reusable sealing unit with flap converter
- choice between single-use dilating trocar pins and reusable standard trocar pins

Cost effective

- an exceptional price-performance ratio leads to a high return on investment

The hybrid trocar solution

Sealing units

Reusable

Single-use, sterile



EKO80P

Sealing unit
for trocars Ø 5 mm
20 units




EK001SU

Single-use sealing unit
for trocars Ø 5 mm
20 units

Trocar sleeves












Reusable



Trocar pins

Reusable

Single-use, sterile



The hybrid trocar solution

Sealing units


Reusable

Single-use, sterile



EKO83P

Sealing unit for trocars
Ø 10/12 mm with
5 mm flap converter



EKO86P

Sealing unit for trocars
Ø 10/12 mm



EKO02SU

Single-use universal
sealing unit for trocars
Ø 10/12 mm

20 units

Trocar sleeves

Reusable



EKO24R

smooth,
with stopcock
110 mm



EKO25R

smooth,
without
stopcock
110 mm



EKO26R

threaded,
with stopcock
110 mm



EKO27R

threaded,
without
stopcock
110 mm

110 mm

110 mm

110 mm

110 mm

150 mm




EKO28R

smooth,
with stopcock
150 mm

Trocar pins

Reusable

Single-use, sterile



The hybrid trocar solution

Sealing units


Reusable

Single-use, sterile



EK083P

Sealing unit for trocars
Ø 10/12 mm with
5 mm flap converter



EK086P

Sealing unit for trocars
Ø 10/12 mm



EK002SU

Single-use universal
sealing unit for trocars
Ø 10/12 mm

20 units

Trocar sleeves

Reusable



EK034R
smooth,
with stopcock
110 mm



EK035R
smooth,
without
stopcock
110 mm



EK036R
threaded,
with stopcock
110 mm



EK037R
threaded,
without
stopcock
110 mm

110 mm

110 mm

110 mm

110 mm


150 mm

EK038R
smooth,
with stopcock
150 mm

Trocar pins

Reusable

Single-use, sterile



The hybrid trocar solution

Trocar for mini-laparotomy (Hasson technique)

Single-use, sterile




EKO02SU

Single-use universal sealing unit for trocars Ø 10/12 mm


20 units

Reusable



EKO86P

Sealing unit for trocars Ø 10/12 mm



EKO24R

smooth, with stopcock 110 mm




EKO64R

blunt (HASSON)



EKO98R

Conical attachment for 10 mm trocar sleeve



All Aesculap endoscopy products on our internet homepage

...just a click away!

Reducing sleeves

Reusable



EK090R

Reducing sleeve,
Ø 10 mm
incl. sealing cap EK093P
5 mm guide



EK091R

Reducing sleeve,
Ø 12 mm
incl. sealing caps
EK093P
5 mm guide
EK094P
10 mm guide



Accessories



EK093P

Sealing cap
5 mm guide
20 units




EK094P

sealing cap
10 mm guide
20 units

Suture anchorage

Reusable



EK097R

Suture fixation for
trocars Ø 5 mm

Sealing units

Reusable



EK080P

Sealing unit for trocars
Ø 5 mm
20 units




consisting of:

Sealing cap EK082P
20 units




Cross slit valve
EK081P
20 units




EK086P

Sealing unit for trocars
Ø 10 - 12 mm
with sealing cap



consisting of:
Sealing cap
EK085P
20 units



Sealing housing
EK088P
1 unit



Cross slit valve
EK084P
20 units



EK083P

Sealing unit for trocars
Ø 10 - 12 mm with
5 mm flap converter




consisting of:
Reducing converter
EK087P
5 units



Sealing housing
EK088P
1 unit



Cross slit valve
EK084P
20 units



Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 74 61 95-0 | Fax +49 (0) 74 61 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

The main product trademark 'Aesculap' is registered trademark of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. C62102

0413/1.0/7